

Tracheal Tube for CAT



Product Features

- Grade soft gel material can provide a safe intubation mechanism, ensuring the comfort and safety of sick animals, and can also accelerate recovery after anesthesia.
- Convenient and safe operation makes the staff more adept when intubating sick animals. The excellent airtight effect can effectively prevent the leakage of anesthetic gas.
- This not only reduces the rejection reaction of sick animals to anesthetic gas but also ensures the life safety of sick animals during anesthesia.
- The rectifier linear design can effectively avoid the occurrence of dead cavities and is more conducive to intraoperative monitoring.
- The anti-bite design can prevent the device from being bitten and blocked by sick animals, preventing gas leakage.

Product size





Dog Tracheal Tube

- A special tracheal tube designed for dog.
- This dog tracheal tube is used to deliver anesthetic gas and/or oxygen to patients.
- Designed to sit in the pharynx, and can avoid many problems caused by tubes such as laryngeal and tracheal trauma.

Dog Tracheal Tube



Feature and Benefits

- Suitable for small dogs and adult large dogs
- Easy to place-particularly useful in emergencies
- No airway narrowing, no increase in airway resistance and protect airway from reflux
- Low dead-space connector
- Excellent seal pressure allowing assisted ventilation but no environmental pollution of anesthetic gases
- No tracheal trauma and no laryngeal damage/spasms
- Safer, smoother recoveries and no post-operative gagging or coughing



Step1: Choose the correct tube.



Step2: Use a laryngoscope to examine the mouth and pharynx, then remove any visible foreign objects or excessive secretions.



Step3: Grasp the tongue with a dry gauze to open the mouth and then insert the tube gently into the pharynx (sterile lube can be used on the edges of tube to slide easily).



Step4: Secure the tie gauze around the tube and tie to pet. Lay the pet down and connect the tube to the anesthesia machine.

| Head shape - Brachycephalic | Lean weight of dog | Head shape - Mesocephalic | Lean weight of dog | Head shape - Dolichocephalic | Lean weight of dog | Size |
|-----------------------------|--------------------|---------------------------|--------------------|------------------------------|--------------------|------|
| | 6-10kg | | 4-7kg | | 2-4kg | D3 |
| | 9-15kg | | 6-10kg | | 4-7kg | D4 |
| | 14-20kg | | 9-15kg | | 6-10kg | D5 |
| | 19-26kg | | 14-20kg | | 9-15kg | D6 |

Remark: this is only for your reference. The appropriate size should be determined by the professional veterinarians based on the practical situation.