



Dog Tracheal Tube

- A special tracheal tube designed for dog.
- This dog tracheal tube is used to deliver anesthetic gas and/or oxygen to patients.
- Designed to sit in the pharynx, and can avoid many problems caused by tubes such as laryngeal and tracheal trauma.

Dog Tracheal Tube



Feature and Benefits

- Suitable for small dogs and adult large dogs
- Easy to place-particularly useful in emergencies
- No airway narrowing, no increase in airway resistance and protect airway from reflux
- Low dead-space connector
- Excellent seal pressure allowing assisted ventilation but no environmental pollution of anesthetic gases
- No tracheal trauma and no laryngeal damage/spasms
- Safer, smoother recoveries and no post-operative gagging or coughing



Step1: Choose the correct tube.






Step2: Use a laryngoscope to examine the mouth and pharynx, then remove any visible foreign objects or excessive secretions.



Step3: Grasp the tongue with a dry gauze to open the mouth and then insert the tube gently into the pharynx (sterile lube can be used on the edges of tube to slide easily).



Step4: Secure the tie gauze around the tube and tie to pet. Lay the pet down and connect the tube to the anesthesia machine.

Head shape - Brachycephalic	Lean weight of dog	Head shape - Mesocephalic	Lean weight of dog	Head shape - Dolichocephalic	Lean weight of dog	Size
	6-10kg		4-7kg		2-4kg	D3
	9-15kg		6-10kg		4-7kg	D4
	14-20kg		9-15kg		6-10kg	D5
	19-26kg		14-20kg		9-15kg	D6

Remark: this is only for your reference. The appropriate size should be determined by the professional veterinarians based on the practical situation.