



## GreenLine®/D™ Fiber Optic Laryngoscopes Disposable Blades & LED Handles



### GreenLine®/D™ Blades

- Constructed of surgical stainless steel
- LED illumination for exceptional, clear viewing
- Soft matte finish prevents reflection and glare during tissue visualization
- Durable nylon blade heel
- Extended safety heel prevents blade-to-handle contact

### GreenLine®/D™ LED Handles

- Fully contained unit
- Textured surface ensures secure grip
- Mercury-free, pre-installed batteries



5-5338-11

5-5338-10

5-5338-09

5-5338-12

### GreenLine®/D™ LED Handles

Item	Size	Description
5-0236-69	Large	Ribbed Grip
5-0236-70	Penlite	Smooth Grip
5-5338-09	Medium	Textured Grip
5-5338-10	Penlite	Textured Grip
5-5338-11	Stubby	Textured Grip
5-5338-12	Standard	Textured Grip



5-0236-69

5-0236-70

GreenLine®/D™ Blades

One time. One choice. One solution.

[www.SunOneSystem.com](http://www.SunOneSystem.com)