


## Close Suction Catheter + Breathing Circuit



### Ultra Unique I Closed Suction Catheter Set

Reduces work of breathing and maintains breathing efficiency

Decreases bacterial load inside of the device

Isolated cleaning peep seal for a cleaner catheter tip

360° Swivel connector with port make more flexible

Ultra clear breathing circuit

Polyurethane sleeve ensures smooth and safe advancement of the catheter



# Closed Suction Catheter Set

## 모델별 제품 규격

모델명	Closed Suction Catheter 규격	멸균여부
L12-5	12Fr	멸균
L14-5	14Fr	멸균

## 구성품 내역

구성품	제품사진	규격	멸균여부
Breathing Circuit		180 cm	멸균
Anesthetic Mask		#5	멸균
Bagging Tube		직경 : 22mm 길이 : 180cm	멸균
HMEF Filter		6(가로) * 6(세로) * 8(높이) (단위: cm)	멸균
Reservoir Bag		용량 : 3L	멸균
CO2 Gas Sampling Line		직경 : 3 mm 길이 : 3 M	멸균